



WRINKLE REMOVAL MACHINE





How Does the 7d hifu fat removal machine? Introduction

High Intensity Focused Ultrasound (HIFU) directly delivers heat energy to skin and subcutaneous tissue that can stimulate and renew the skin's collagen and thus consequently improving the texture and reducing sagging of the skin. It literally achieves the results of a facelift or a body lifts without any invasive surgery or injections, moreover, an added bonus of this procedure is that there is no downtime.



How HIFU Skin Tightening works?

Introduction

7D technology has a smaller focus point than other. With more accurately transmits high-energy focused ultrasound energy at 65~75°C to the target skin tissue layer, which results in a hot coagulation effect without harming surrounding tissues. While stimulating the proliferation of collagen and elastic fibers, it greatly improves the comfort and give you a perfect face with the skin firm, and face tightening


HIFU(High Intensity Focused Ultrasound) is delivered into the deep skin fractionally



Cause thermal coagulation to SMAS (Superficial Muscular Aponeurotic System), fascia, deep dermis and subcutaneous fibrous tissue.




By wound healing, skin tightening and rejuvenation happens gradually with time



Whole body Treatment:





Seven kinds of Depth treatment probe

Macro focused transducers with a larger density and deeper penetration of the technology target the body's subcutaneous layers at 3 different depths to slim contours and remove stubborn fat. Multiple cartridges are customized for shallow layers of the face and sagging skin, collagen renewal, and wrinkle reduction.



01

Lift & Tighten


Four multi-depth cartridges directly transfer ultrasound energy into multiple layers of the face to trigger collagen renewal for a tighter , younger complexion . The 2.0 mm cartridge can also treat both narrow and sensitive contours with a slimmer transducer and faster shot speed



Tighten & Contouring body


Macro-focused cartridges at 6.0mm, 9.0mm, and 13.0mm directly transfer an equal distribution of ultrasound energy into the subcutaneous tissue of multiple body areas, inducing contractions as well as naturally disposing of stubborn fat cells for contouring purposes








Lasting And More Effective



Real picture show









7d hifu technology wrinkle removal machine

face liftin+skin tightening+body contouring

7D ULTRAMAGE



System Specification V	
7D	
7d hifu technology wrinkle	
Removal machine	
220V 50HZ	
150W	
2-7MHz	
1.0mm 2.0mm 3.0mm,4.5mm	
6.0mm,9.0mm,13mm	
white	